

# LegholderRX® C €

The „Must“ for Surgery  
of the Knee Joint!



# LegholderRX® C €



Used in almost every hospital in Switzerland and in more than 30 countries throughout the world.

Made in Switzerland, Swissmedic CH-201505-0035; Registered with the FDA

**Can be used on all common  
operation tables**





# Delivered with Velcro belts



# 6 different stable positions



**In combination with a lateral post at the middle of the thigh one assistant become overfluous**





# X-ray transparent




**Position: The highest support just proximal of the insertion of the achilles tendon**





# Lateral post at the middle of the thigh



# 45° Flexion



# 80° Flexion





# 110° Flexion



# 130° Flexion



# 150° Flexion






**After final installation →  
usual disinfection**



**steril prepping**



# finished prepping





**Replace the lateral post under the prepping at the middle of the thigh**



# 45° Flexion



# 80° Flexion





# 110° Flexion



# 130° Flexion



# 150° Flexion





# LegholderRX® CE



**The „Must“ for Surgery of the Knee Joint!**

Wirz Ortho Mechanics

Bernstrasse 1

CH-3076 Worb

Tel: +41318324353

Fax +41318324354

Mobile +41796306644

[mechanik@wirz-ortho.ch](mailto:mechanik@wirz-ortho.ch)

[www.wirz-ortho.ch](http://www.wirz-ortho.ch)

Swissmedic CH-201505-0035

Registered with the FDA

