


Rolimeter[®] € P/N 50A Developed by Roland Jakob, Motier, CH

The Rolimeter is designed for measuring anterior and posterior knee joint laxity during Lachman, anterior drawer, and “step off” tests



Lachman testing

- The second hand pushes the feeler together with the tibia down
- The calf is pulled up with the second hand
- The drawer can be measured



**WIRZ ORTHO
MECHANICS**

MEDIZINTECHNIK UND FEINMECHANIK
Bernstrasse 1, CH-3076 Worb, Tel: +41 31 832 43 53
Werkstatt: Hämli mattweid 101, CH-3508 Ami,
Mobile: Peter Wirz +41 79 630 66 44, Fax: +41 31 832 43 54
mechanik@wirz-ortho.ch www.wirz-ortho.ch